

## ALLFIT™

### A. GENERAL PRODUCT INFORMATION


ALLFIT are low temperature thermoplastic pellets that can be kneaded and formed to any shape. This product is used in physical rehabilitation to make individual web spacers, add-ons to orthoses, to fabricate customized handles, utensils, pencils and all kinds of tools. The pellets can also be used for desensitization exercises.

! ALLFIT pellets are not suitable for internal use. They may not be used on open wounds or in the mouth.


### B. PRODUCT RANGE

Art. No.	Product description
703011	Packaging: 6 bags of 500g

### C. ACTIVATION AND APPLICATION TECHNIQUE




1. Pour ALLFIT in a Suspan (Art.no. 35240), a water heater or a bowl of hot water at 65°C (149°F) or more.



2. The pellets will stick together and become transparent to form a viscous mass.



3. Knead them homogeneously.



4. Form it around an object into the desired shape. ALLFIT sticks to most orthotic fabrication materials that are not coated and to porous surfaces.
5. Cooling can be enhanced by means of cold water or cold spray.
6. ALLFIT is completely or partially remouldable by reheating the moulded ALLFIT in hot water of 65°C (149°F) it is possible to correct and reshape the ALLFIT into the desired form or to remove it from the chosen object and apply onto another object.

**Never expose an object adapted with ALLFIT to temperatures higher than 60°C (140°F) ( e.g. behind the windshield of a car in the sun, in the proximity of a hot radiator or a fireplace).**

#### D. PRECAUTIONS BEFORE USE

High temperatures up to a maximum of 120°C (248°F) do not damage ALLFIT pellets, but are not user-friendly.

#### E. MAINTENANCE AND WASTE MANAGEMENT

ALLFIT can be disposed of with normal household waste without harming the environment. ALLFIT is biodegradable.

#### F. STORAGE

After use, ALLFIT must be stored in its original packaging and in a dry place at a temperature not lower than 10°C (50°F) or higher than 30°C (86°F).

Low temperature thermoplastics can only be kept for a limited period of time and must be protected as much as possible from light, heat and humidity. The material ages in direct relation to storage circumstances.

#### G. GENERAL SAFETY ADVICE

Check the temperature of mouldable ALLFIT before adapting it to the patient's hand.

**!• ALLFIT is not suitable for internal use and may not be used on open wounds or in the mouth.**

**!• Never use an open flame to activate ALLFIT pellets.**

#### H. ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets and regulatory information, please visit our website [www.orfit.com](http://www.orfit.com).

#### Note:

The instructions were written in accordance with the European Directive 93/42/EEC for Medical Devices. It is prohibited to make alterations to this text without prior approval from Orfit Industries. ALLFIT™ is a trademark from ORFIT INDUSTRIES N.V.



**ORFIT INDUSTRIES**  
Vosveld 9A | B-2110 Wijnegem | Belgium  
T (+32) (0)3 326 20 26  
[welcome@orfit.com](mailto:welcome@orfit.com)

[www.orfit.com](http://www.orfit.com)



Ref. No. 31150  
VERSION 3  
LAST UPDATE: 25/01/2019  
REVISION DATE: 25/01/2021