

PUSH PIN UNLOCK TOOL



Article No. : 32199

A. GENERAL PRODUCT INFORMATION

The product referred to in these instructions is a tool that can be used to unlock push pins and remove them from a push pin base plate used for patient positioning and immobilization in radiation therapy.

B. PRODUCT DESCRIPTION

This Push Pin Unlock Tool is specifically designed to easily unlock and release push pins from a push pin base plate.

C. HOW TO USE?

To unlock the push pin, first slide the tool under the pin as shown in the pictures below. Place your finger at the opposite end of the tool and press down until the pin is lifted.



D. STORAGE

Always store the product in a safe place to prevent it from getting damaged.

Store the product between +10°C (50°F) and 30 °C (86°F).

E. PROPERTIES

The following physical properties apply to this product:

Dimensions: L 8 cm x W 2.5 cm x H 1 cm
L 3.15" x W 0.98" x H 0.39"

Weight: 20 g
0.04 lbs

Material : Aluminum

F. MAINTENANCE AND WASTE MANAGEMENT

This product can be cleaned and disinfected by means of soapy water or an isopropanol based disinfectant, applied with a soft cloth. If unsure about the cleaning fluid, do not use. **Never use aerosol sprays, corrosive cleaning agents, solvents or abrasive detergents.**

Periodic checks of the product should be done to insure the parts are not worn and require repair or replacement. **Do not attempt to make repairs yourself.** Contact your distributor if there are any questions or concerns.

The product can be disposed of with household waste.

G. ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets and regulatory information, please visit our website www.orfit.com.

Note:

The instructions were written in accordance with the European Directive 93/42/EEC for Medical Devices. It is prohibited to make alterations to this text without prior approval from Orfit Industries.



ORFIT INDUSTRIES
Vosveld 9A | B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcome@orfit.com

www.orfit.com



Ref. No. 50185
VERSION 1
LAST UPDATE: 10/01/2019
REVISION DATE: 10/01/2021