


EFFICAST® HIGH PRECISION MASKS FOR HEAD, NECK AND SHOULDERS

4-POINT MASK - EFFICAST®



Art. No. 35764/2MA Thickness: 2 mm. - Perforation: maxi.	Art. No. 35764/16MI Thickness: 1,6 mm. - Perforation: micro.
--	--

Application: Head & Shoulders supine.
Material: Efficast®.
Type: 4-point Head, Neck & Shoulders Mask with nose hole.
Packaging: 10 per box.



Art. No. 35764/2MA/NH Thickness: 2 mm. - Perforation: maxi.

Application: Head & shoulders supine.
Material: Efficast®.
Type: 4-point Head, Neck & Shoulders Mask without nose hole.
Packaging: 10 per box.



Art. No. 35764/2MA/EM Thickness: 2 mm. - Perforation: maxi.	Art. No. 35764/16MI/EM Thickness: 1,6 mm. - Perforation: micro.
---	---

Application: Head & Shoulders supine. For patients who feel claustrophobic and for application of a bite block or canula.
Material: Efficast®.
Type: 4-point Head, Neck & Shoulders Mask with nose, eye and mouth hole.
Packaging: 10 per box.



Art. No. 35764/2MA/M Thickness: 2 mm. - Perforation: maxi.	Art. No. 35764/16MI/M Thickness: 1,6 mm. - Perforation: micro.
--	--

Application: Head & Shoulders supine. For application of a bite block or canula.
Material: Efficast®.
Type: 4-point Head, Neck & Shoulders Mask with nose and mouth hole.
Packaging: 10 per box.

5-POINT MASK - EFFICAST®



Art. No. 33700/2MA Thickness: 2 mm. - Perforation: maxi.	Art. No. 33700/16MI Thickness: 1,6 mm. - Perforation: micro.
--	--

Application: Head & Shoulders supine. For IMRT.
Material: Efficast®.
Type: 5-point Head, Neck & Shoulders Mask with nose hole.
Packaging: 10 per box.



Art. No. 33700/2MA/EM Thickness: 2 mm. - Perforation: maxi.	Art. No. 33700/16MI/EM Thickness: 1,6 mm. - Perforation: micro.
---	---

Application: Head & Shoulders supine. For IMRT. For patients who feel claustrophobic and for application of a bite block or canula.
Material: Efficast®.
Type: 5-point Head, Neck & Shoulders Mask with nose, eye and mouth hole.
Packaging: 10 per box.



Art. No. 33700/2MA/M Thickness: 2 mm. - Perforation: maxi.	Art. No. 33700/16MI/M Thickness: 1,6 mm. - Perforation: micro.
--	--

Application: Head & Shoulders supine. For IMRT. For application of a bite block or canula.
Material: Efficast®.
Type: 5-point Head, Neck & Shoulders Mask with nose and mouth hole.
Packaging: 10 per box.



Art. No. 33729/2MA Thickness: 2 mm. - Perforation: maxi.
--

Application: Head and Shoulders supine. For IMRT.
Material: Efficast®.
Type: 5-point Head, Neck & Shoulders. For large patients that need extra material around the shoulders and over the chest.
Packaging: 10 per box.



Art. No. 33730 Thickness: 2 mm. - Perforation: maxi.
--

Application: Head & Shoulders Supine for ultra-high precision treatments (SRT, SRS, Hypofractionation, Single dose Treatments)
Material: Efficast® 2 mm maxi, reinforced with low density Orflight™ 1,6 mm micro.
Packaging: 10 per box.