


PRE-CUTS FOR HEAD, NECK AND SHOULDER FIXATION

ON OTHER HARDWARE

Now you can realize the superior patient immobilization performance of Efficast® when used with your fixation hardware.

Choose from a wide variety of Efficast® pre-cut masks with attached profiles that have been specifically designed to provide superior moulding and fixation properties and superior patient comfort.



Art. No. 33690/2MA
Pre-cut with profiles

Art. No. 33690/2MA/NP
Pre-cut without profiles

Application: Head supine
Material: Efficast®
Type: 3-point pediatric mask
Thickness: 2 mm
Perforation: maxi
Packaging: 10 per box



Art. No. 33770/2MA
Pre-cut with profiles

Art. No. 33770/2MA/NP
Pre-cut without profiles


Application: Head, neck & shoulders supine - IMRT
Material: Efficast®
Type: 5-point mask for regular sized patients
Thickness: 2 mm
Perforation: maxi
Packaging: 10 per box



Art. No. 33770/2MI
Pre-cut with profiles

Art. No. 33770/2MI/NP
Pre-cut without profiles

Application: Head, neck & shoulders supine - IMRT
Material: Efficast®
Type: 5-point mask for regular sized patients
Thickness: 2 mm
Perforation: micro
Packaging: 10 per box



Art. No. 35772/2MA
Pre-cut with profiles

Art. No. 35772/2MA/NP
Pre-cut without profiles

Application: Head supine
Material: Efficast®
Type: 2-point mask
Thickness: 2 mm
Perforation: maxi
Packaging: 10 per box



Art. No. 35772/16MI
Pre-cut with profiles

Art. No. 35772/16MI/NP
Pre-cut without profiles

Application: Head supine
Material: Efficast®
Type: 2-point mask
Thickness: 1.6 mm
Perforation: micro
Packaging: 10 per box



Art. No. 35773/2MA
Pre-cut with profiles

Art. No. 35773/2MA/NP
Pre-cut without profiles

Application: Head supine
Material: Efficast®
Type: 3-point mask
Thickness: 2 mm
Perforation: maxi
Packaging: 10 per box