



Art. No. 35773/16MI
Pre-cut with profiles

Art. No. 35773/16MI/NP
Pre-cut without profiles

Application: Head supine
Material: Efficast®
Type: 3-point mask
Thickness: 1.6 mm
Perforation: micro
Packaging: 10 per box



Art. No. 35773/2MI
Pre-cut with profiles

Art. No. 35773/2MI/NP
Pre-cut without profiles

Application: Head supine
Material: Efficast®
Type: 3-point mask
Thickness: 2 mm
Perforation: micro
Packaging: 10 per box



Art. No. 35774/2MA
Pre-cut with profiles

Art. No. 35774/2MA/NP
Pre-cut without profiles

Application: Head, neck & shoulders supine
Material: Efficast®
Type: 4-point mask
Thickness: 2 mm
Perforation: maxi
Packaging: 10 per box



Art. No. 35774/16MI
Pre-cut with profiles

Art. No. 35774/16MI/NP
Pre-cut without profiles

Application: Head, neck & shoulders supine
Material: Efficast®
Type: 4-point mask
Thickness: 1.6 mm
Perforation: micro
Packaging: 10 per box



Art. No. 35791/2MA
Pre-cut with profiles

Art. No. 35791/2MA/NP
Pre-cut without profiles


Application: Head, neck & shoulders supine - IMRT
Material: Efficast®
Type: 5-point mask for large sized patients
Thickness: 2 mm
Perforation: maxi
Packaging: 10 per box



Art. No. 35791/2MI
Pre-cut with profiles

Art. No. 35791/2MI/NP
Pre-cut without profiles

Application: Head, neck & shoulders supine - IMRT
Material: Efficast®
Type: 5-point mask for large sized patients
Thickness: 2 mm
Perforation: micro
Packaging: 10 per box



Art. No. 35791/16MI
Pre-cut with profiles

Art. No. 35791/16MI/NP
Pre-cut without profiles

Application: Head, neck & shoulders supine - IMRT
Material: Efficast®
Type: 5-point mask for large sized patients
Thickness: 1.6 mm
Perforation: micro
Packaging: 10 per box