


# LOW DENSITY HEAD SUPPORTS

High Precision head supports with a 3-dimensional shape for comfortable, stable and exact positioning of the head. The cranial stop precisely cradles the skull, resulting in an excellent reference point for a highly reproducible position.


Made of low density coated foam with a low attenuation factor.  
Easy to clean.



Art. No. 35755  
Size A with neck flaps  
Fits to 90% of all patients.




Art. No. 35765  
Size A without neck flaps  
Fits to 90% of all patients.



Art. No. 35753  
Size B with neck flaps 80%  
scaled down version of size A,  
for smaller sized patients.



Art. No. 35753ZF  
Size B without neck flaps  
80% scaled down version  
of size A, for smaller sized  
patients.



Art. No. 35752  
Size C with neck flaps  
For patients with a short neck.



Art. No. 35752ZF  
Size C without neck flaps  
For patients with a short  
neck.



Art. Nr. 32704  
Set of the above 6 low density  
head supports.

