


# HIGH PRECISION HARDWARE SYSTEM FOR HEAD, NECK AND SHOULDERS



Art. No. 32110  
High Precision Base plate.

Application: head, neck and shoulders fixation.  
Material: Carbon Fibre sandwich construction with Rohacel™ foam.  
Width: 54 cm.



Art. No. 32150  
High Precision Long Base plate for IMRT.

Application: head, neck and shoulders.  
Material: Carbon Fibre sandwich construction with Rohacel™ foam.  
Width: 54 cm.




Art. No. 32111  
High Precision Pediatric Base Plate.

Application: head and Neck.  
Material: Carbon Fiber sandwich construction with Rohacell® foam.  
Width: 34 cm



Art. No. 32130  
High Precision Small Base Plate for IMRT.

Application: head only.  
Material: Carbon Fibre sandwich construction with Rohacel™ foam.  
Width: 19.9 cm.



Art. No. 32140  
High Precision Secondary Base Plate.

Application: head only when used with 40 mm block and 18° wedge.  
Material: Carbon Fibre sandwich construction with Rohacel™ foam.  
Width: 19.9 cm.



Art. No. 35703  
High Precision block – 20 mm.

Application: elevation of the head.  
Material: Carbon Fibre sandwich construction with Rohacel™ foam.



Art. No. 35706  
High Precision block – 40 mm.

Application: elevation of the head.  
Material: Carbon Fibre sandwich construction with Rohacel™ foam.



Art. No. 35704  
High Precision wedge – 9° angle.

Application: flexion and extension of the head.  
Material: Carbon Fibre sandwich construction with Rohacel™ foam.



Art. No. 35707  
High Precision wedge – 18° angle.

Application: flexion and extension of the head.  
Material: Carbon Fibre sandwich construction with Rohacel™ foam.



Art. No. 32700  
Set of art. No. 35703, 35704, 35706, 35707.

Material: Carbon Fibre sandwich construction with Rohacel™ foam.



Art. No. 32110-COMBI  
Set of art. No. 32110, 35703, 35704, 35706, 35707.

Material: Carbon Fibre sandwich construction with Rohacel™ foam.



Art. No. 32150-COMBI  
Set of art. No. 32150, 35703, 35704, 35706, 35707.

Material: Carbon Fibre sandwich construction with Rohacel™ foam.



Art.No. 32380  
Prone Head Support (adjustable).  
Application: prone position of the head.

Material: Low density foam- coated. Adjustable system on high density plate.