

UON™ - DUON™ HEAD AND SHOULDER FIXATION

This two-part system will eliminate the inconvenience and difficulty you experience when placing and removing the very large first generation S-Shaped mask from the water bath. Moulding DUON™ to your patients will be quicker, easier and more convenient. DUON™ eliminates the need for you to purchase and inventory separate U-Shaped head masks and separate full size S-Shaped Masks saving you money.

This new and innovative approach combining a separate U-Shaped and S-Shaped mask in one product has been pioneered by Orfit Industries n.v.

DUON™ is designed to be compatible for use with your S-Shaped hardware.

Orfit Industries also offers its own ultra light weight carbon fibre laminate and high pressure laminate base plates for use with DUON™.



Art. No. 33101
Integra Type-S board for Head and neck immobilization with UON and DUON and with Efficast HP Head and neck masks

DUON™ is a U-Shaped mask designed for superior head patient fixation. But DUON™ is more. It is also an S-Shaped mask for superior neck and shoulder fixation. When interlocked together using the quick lock push pins the two DUON™ parts instantaneously become a superior S-Shaped head and shoulder patient fixation mask. Moulding superior head, head and shoulder patient fixation masks has never been quicker, easier and more cost-effective.



DUON Head only Baseplate in Low density Carbon Fiber material.
Art. No. 32095



DUON Head Only Baseplate in High density material
Art. No. 32090



Art. No. 60095
Thickness: 3.2 mm.
Perforation: Maxi
Application: Head only
Material: U-Plast High Stretch
Packaging: 10 per box



Art. No. 60097
Thickness: 3.2 mm.
Perforation: Maxi
Application: Head only with additional material over the chin
Material: U-Plast High Stretch
Packaging: 10 per box



Art. No. 60117
Thickness: 2.4 mm.
Perforation: Micro Plus
Application: Head only with additional material over the chin
Material: U-Plast High Stretch
Packaging: 10 per box



Art. No. 60120
Thickness: 2.4 mm.
Perforation: Hybrid
Application: Head only
Material: U-Plast High Stretch
Packaging: 1 per box



Art. No. 60122
Thickness: 2.4 mm.
Perforation: Micro Plus
Application: Head only with additional material over the chin
Material: U-Plast High Stretch
Packaging: 1 per box



Art. No. 60110
Thickness: 2.4 mm.
Perforation: Hybrid
Application: Head only
Material: U-Plast High Stretch
Packaging: 10 per box



Art. No. 60112
Thickness: 2.4 mm.
Perforation: Micro Plus
Application: Head only
Material: U-Plast High Stretch
Packaging: 10 per box



Art. No. 60130
Thickness: 2.0 mm.
Perforation: Maxi
Application: Shoulder immobilization in combination with DUON Head
Material: U-Plast Low Stretch
Size: Small
Packaging: 1 per box



Art. No. 60140
Thickness: 2.4 mm.
Perforation: Maxi
Application: Shoulder immobilization in combination with DUON Head
Material: U-Plast Low Stretch
Size: Large
Packaging: 1 per box



Art. No. 60123
Thickness: 3.2 mm.
Perforation: Maxi
Application: Head only with additional material over the chin
Material: U-Plast High Stretch
Packaging: 1 per box