

PELVICAST™ HIGH PRECISION BELLY AND ABDOMEN PATIENT FIXATION SYSTEM



Art. No. 32122
High Precision Base plate.
Application: Pelvic fixation with 4-points masks
Material: Carbon Fiber sandwich construction with Rohacel foam
Width: 54 cm



Art. No. 32123
High Precision Base plate.
Application: Pelvic fixation with 4-points masks
Material: High density High Pressure laminate
Width: 54 cm



Art. No. 35781-N1/2
Central panel for added patient support and comfort.
Application: supine position.
Material: High Pressure Laminate sandwich construction with Rohacel™ foam.



Art. No. 35782
Central panel for added patient support and comfort.
Application: supine position.
Material: Carbon Fibre Laminate.




Art. No. 32037
Leg separator for 2nd generation Pelvicast™ and AIO Base plates.
Contains fixation slots for 6-point Pelvicast™ masks.
Material: POM.



Art. No. 32122-COMBI1
High precision base plate with Leg separator and fixator
Application: Pelvic fixation with 6-points masks
Material: Carbon Fiber sandwich construction with Rohacel foam
Width: 54 cm

Art. No. 32123-COMBI1
High precision base plate with Leg separator and fixator
Application: Pelvic fixation with 6-points masks
Material: High density High Pressure laminate.
Width: 54 cm



Art. No. 32122-COMBI2
High precision base plate with Leg separator and fixator
Small base plate with head support and Tablefix Leg and Knee Cushion Set (THE AIO SOLUTION®)
Application: Very high precision Pelvic fixation with head and knee indexing.
Material: Carbon Fiber pelvic base plate with high density head base plate and head support.

Art. No. 32123-COMBI2
High precision base plate with Leg separator and fixator
Small base plate with head support and Tablefix Leg and Knee Cushion Set (THE AIO SOLUTION®)
Application: Very high precision Pelvic fixation with head and knee indexing.
Material: High density pelvic base plate with high density head base plate and head support.

