

# PELVICAST™ HIGH PRECISION BELLY AND ABDOMEN PATIENT FIXATION SYSTEM



Art. No. 35710/32MA  
Application: Pelvis, prone & supine.  
Material: Efficast®.  
Type: 6-point mask.  
Size: Small.  
Thickness: 3.2 mm.  
Perforation: maxi.  
Packaging: 10 per box.



Art. No. 35784/32MA  
Application: Pelvis, prone + supine.  
Material: Efficast®.  
Type: 4-point Pelvis Mask Large, Prone & Supine.  
Size: Large.  
Thickness: 3.2 mm.  
Perforation: maxi.  
Packaging: 10 per box.



Art. No. 35711/32MA  
Application: Pelvis, prone & supine.  
Material: Efficast®.  
Type: 6-point mask.  
Size: Medium.  
Thickness: 3.2 mm.  
Perforation: maxi.  
Packaging: 10 per box.



Art. No. 35780/32MA  
Application: Pelvis, supine.  
Material: Efficast®.  
Type: 6-point Pelvis Mask Large, Supine.  
Size: Large.  
Thickness: 3.2 mm.  
Perforation: maxi.  
Packaging: 10 per box.



Art. No. 35712/32MA  
Application: Pelvis, prone & supine.  
Material: Efficast®.  
Type: 6-point mask.  
Size: Large.  
Thickness: 3.2 mm.  
Perforation: maxi.  
Packaging: 10 per box.



Art. No. 35790/2MA  
Application: Breast supine.  
Material: Efficast®.  
Type: 2-point mask.  
Thickness: 2 mm.  
Perforation: maxi.  
Packaging: 10 per box.

