

BASE PLATE FIXATION DEVICES



Art. No. 32001 - Loxon™ Set for HP System (*High Precision*)
Suitable only for the HP System and for Table Tops with Loxon™ slots.



Art. No. 32002 - Loxon™ Set for Pelvicast™ System
Suitable only for the Pelvicast™ Base Plates and for Table Tops with Loxon™ slots.



Art. No. 35781-N1/6
Pair of marking knobs for indexing a vacuum or other mattress on a base plate.
Material: POM.



Art. No. 35745
Table Fix for all other couch types.
Material: light weight aluminium.



Art. No. 32160
Table Fix for Pelvicast™ System HPL Base Plate.
Set of 4 pcs.
Material: light weight aluminium.



Art. No. 32003
Loxon™ for THE AIO SOLUTION® Base plate.
Suitable only for indexing the AIO 2nd generation Base plate to an Exact® couch.