

AIO HEAD & NECK APPLICATIONS

The main goal in the positioning of patients with head & neck tumors is a very precise and reproducible positioning and immobilization of the head and upper torso of the patient. To achieve this, a combination of the AIO hardware and Efficast® thermoplastic masks provides an ideal solution.

- The AIO base plate allows for a reproducible and comfortable positioning of the patient, reducing stress and anxiety.
- Efficast® 3-points masks assure a precise immobilization of the head and 5-points masks of the head and shoulders.



AIO Head & Neck Applications

