

AIO STORAGE CABINET

The AIO Storage Cabinet has 5 vertical sliding panels, and can store a complete AIO cushion set, as well as all fixation devices and an AIO board.



The cabinet is made of laminated pressed wood boards, which are easy to clean. It is easy to assemble, but at the same time sturdy to make it long-lasting. The cushions are stored in such a way that they remain well protected from possible damage.

The panels slide out completely, allowing easy access to the cushions.

The cabinet is designed to efficiently store the AIO Solution® components in the simulation or treatment room while taking up limited space (width x depth x height = 1636x780x1826 mm)

Every AIO part has its own dedicated place so that it can be easily found. This reduces patient set-up time and makes inventory control easy in the event of misplaced parts.

