



SRS-Fix enhances patient comfort and safety and adds to the value of stereotactic radiotherapy.

SRS-Fix Immobilization

Value for User

Value for Patient

Easy, fast, trouble-free mask preparation and setup

Short mask molding procedure thanks to Orfit's proven single shell multi-point mask

Rapid procedure for optimal patient throughput

Short time spent for simulation

Short daily setup time thanks to Orfit's proven single shell multi-point mask

Increased patient throughput on treatment machine

Short stay on treatment table

Easy-to-apply BiteFix

Reduces pitch movement

Minimal invasive

High technology thermoplastic mask material

Hybrid mask based on nanotechnology with low shrinkage and high fixation force

Excellent immobilization

Increased level of comfort

Antibacterial and non-stick surface

No risk of cross-contamination in the department

Less exposure to bacteria on mask surface and no sticking of the mask to hair and face

High precision and reproducibility

Personalized head support

High precision of positioning

High level of comfort

Mask with high definition of the shape of the head

High precision of positioning

Better clinical effect


System based on existing positioning hardware

Orfit High Precision Base Plate and Head Supports form the basis of the system

Low or inexistant investment cost

Product highlights

A mask that makes a difference: the Orfit SRS-Fix mask



Its hybrid configuration with a Nanor® layer, its controlled shrinkage behavior and strong fixation force result in optimal immobilization with a high degree of patient comfort.

BiteFix

BiteFix offers a simple, one-size-fits-all solution to reduce pitch movement. The airway and mask opening ensures comfortable breathing.



Customized head support and vacuum bag

Create a personalized patient immobilization solution with your preferred choice of individual head rest: Thermofit or Vacuum Bag.



SRS-Fix

SRS treatment puts high demands on patient immobilization. SRS-Fix is simple in use, reliable in functioning and stable in immobilizing a patient's head.

It is a patient-centered solution that fulfills the special requirements to achieve millimeter precision and daily reproducibility.



About ORFIT

Orfit brings high precision and comfort to the positioning and immobilization of cancer patients in Radiation Oncology. The combination of Orfit thermoplastic masks and positioning devices provides optimal stability.

Discover the complete family of Orfit systems for brain, head and neck, supine and prone breast treatment, pelvis and abdomen treatment, SRS, SBRT, extremities, pediatric, MR and proton therapy on www.orfit.com



SRS-Fix Solution

Simple, Reliable and Stable Fixation
for stereotactic head and neck treatment

Visit
WWW.ORFIT.COM
for more information.

www.orfit.com



ORFIT INDUSTRIES
Vosveld 9A
B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcome@orfit.com

ORFIT INDUSTRIES AMERICA
810 Ford Drive
Norfolk, VA 23523 | USA
T 516-935-8500
welcome@orfit.com

01/12/2019
Art. N° 51025

orfit The
precision
expert