



The  
reliable  
choice

IN PROSTHETIC AND ORTHOTIC  
THERMOPLASTIC MATERIALS

— MORE ORFITRANS —



— MORE POSSIBILITIES —

# More Orfitrans More possibilities

Orfitrans is an advanced range of high temperature thermoplastic sheet materials for the production of check sockets, flexible sockets and orthoses.

The superb esthetic qualities and outstanding handling and performance characteristics of Orfitrans will help to eliminate or minimize many of the difficulties that you may experience with the materials you currently use.


The Orfitrans product family consists of different thermoplastic materials that have been specifically developed to provide solutions for a wide range of orthotic and prosthetic applications. Select the Orfitrans material with the performance characteristics that will best meet your specific requirements and enable you to achieve the desired patient outcome.

# ORFITRANS STIFF

---

Orfitrans Stiff is the Gold Standard for test sockets. It will not break, crack or shatter when dropped. The material is virtually indestructible, but you can still easily drill into it.

No other material comes close to Orfitrans Stiff in performance and quality. It pulls easily and uniformly, and the material feels great. Once you start to use Orfitrans Stiff you will immediately know why it is the "GOLD STANDARD".



# ORFITRANS MEDIUM SOFT

---

Orfitrans Medium Soft combines enough rigidity to give control and enough flexibility to allow range of motion. It is structurally sound and will not cold flow at thinner areas of the flexible interface.

The material will maintain its shape and the volume does not change or stretch outward. It allows for the change or flattening of the limb when sitting. This ensures that the socket will conform to the change to prevent loss of suspension while maintaining comfort. It is ideal to use for vacuum transtibial and transfemoral sockets.



Medium Soft in Natural Color

Medium Soft Beige **NEW**

# ORFITRANS EXCEL

---

Orfitrans Excel is lighter in weight than most other materials for flexible prosthetic sockets, resulting in increased patient comfort.

This low friction, flexible material is not tacky nor sticky to mold. The material pulls easily and uniformly. The surface finish is silky and smooth and buffs easily. Both bubble and drape forming are possible. The material will bond to itself permanently.



Excel

Excel Black

Orfitrans Excel will not cold flow and is not subject to discoloration. It is resistant to wrinkling and very stable over time. The moderate bending modulus of Orfitrans Excel allows for functional mobility and muscle movement.

# ORFITRANS EXTRA SOFT

# ORFITRANS EXTRA SOFT SILICONE

---

## ORFITRANS EXTRA SOFT

Orfitrans Extra Soft is recommended for highly flexible prosthetic sockets for active patients. The material is extremely smooth and is used for its softness, patient friendliness and its comfort. Orfitrans Extra Soft is semi-transparent.

Orfitrans Extra Soft Beige comes in an extra thin sheet for the use of upper extremity sockets.



Extra Soft

Extra Soft Beige **NEW**

Extra Soft Silicone

## ORFITRANS EXTRA SOFT SILICONE

Orfitrans Extra Soft Silicone has the same excellent properties as Orfitrans Extra Soft. In addition it has the most silky surface of all Orfitrans products. The product comes in a non-transparent color.

# ORFITRANS SUPRA SOFT **NEW**

---

Orfitrans Supra Soft is an extremely flexible, silky-smooth non-stick material. It is the softest sheet material available in the Orfitrans product range and the product's optimum comfort makes it ideal for flexible prosthetic sockets for geriatric patients.



Supra Soft **NEW**

---

Find fresh insights, education and updates of expertise on the Orfit O&P blog  
[www.orfit.com/prosthetics-orthotics/blog](http://www.orfit.com/prosthetics-orthotics/blog)

Go to [www.orfit.com](http://www.orfit.com) to consult article numbers,  
instructions for use and other product information.



[www.orfit.com](http://www.orfit.com)



**ORFIT INDUSTRIES**

Vosveld 9A | B-2110 Wijnegem | Belgium  
T (+32) (0)3 326 20 26  
[welcome@orfit.com](mailto:welcome@orfit.com)

**ORFIT® INDUSTRIES AMERICA**

810 Ford Drive | Norfolk, Virginia 23523  
T 516-935-8500  
[welcome@orfit.com](mailto:welcome@orfit.com)